


Every Day. Every Case.



ARTIVION™

E-tegra™
Stent Graft System

Ordering Information

E-tegra Stent Graft System

01 Main Body



Catalog Number Bifurcation Prosthesis	Proximal Ø (mm)	Distal Ø (mm)	Total Trunk Length (mm)	Covered Length (mm)	OD Delivery System (F/mm)
93MB2313L10-08	23	13	85	100	18/6.1
93MB2313L13-10	23	13	100	125	18/6.1
93MB2613L10-08	26	13	85	100	18/6.1
93MB2613L13-10	26	13	100	125	18/6.1
93MB2913L10-08	29	13	85	100	18/6.1
93MB2913L13-10	29	13	100	125	18/6.1
93MB3213L10-08	32	13	85	100	20/6.8
93MB3213L13-10	32	13	100	125	20/6.8
93MB3613L10-08	36	13	85	100	20/6.8
93MB3613L13-10	36	13	100	125	20/6.8

Special sizes on demand only

93MB2313L15-08	23	13	85	155	18/6.1
93MB2613L13-08	26	13	85	135	18/6.1
93MB2616L15-08	26	16	85	155	18/6.1
93MB2913L15-08	29	13	85	155	18/6.1
93MB2919L15-10	29	19	100	160	18/6.1
93MB3216L15-08	32	16	85	155	20/6.8
93MB3219L15-10	32	19	100	160	20/6.8

02 Contralateral Leg



Catalog Number	Proximal Ø (mm)	Distal Ø (mm)	Usable Length (mm)	Total Length (mm)	OD Delivery System (F/mm)
93CL1510L05	15	10	50	80	16/5.7
93CL1510L07	15	10	65	95	16/5.7
93CL1510L09	15	10	90	120	16/5.7
93CL1510L10	15	10	100	130	16/5.7
93CL1513L05	15	13	50	80	16/5.7
93CL1513L07	15	13	70	100	16/5.7
93CL1513L09	15	13	90	120	16/5.7
93CL1513L10	15	13	105	135	16/5.7
93CL1516L03	15	16	30	60	16/5.7
93CL1516L05	15	16	50	80	16/5.7
93CL1516L07	15	16	70	100	16/5.7
93CL1516L09	15	16	90	120	16/5.7
93CL1516L10	15	16	105	135	16/5.7
93CL1519L05	15	19	50	80	16/5.7
93CL1519L07	15	19	70	100	16/5.7
93CL1519L09	15	19	90	120	16/5.7
93CL1519L10	15	19	105	135	16/5.7
93CL1522L05	15	22	50	80	16/5.7
93CL1522L07	15	22	70	100	16/5.7
93CL1522L09	15	22	90	120	16/5.7
93CL1522L10	15	22	105	135	16/5.7
93CL1525L05	15	25	50	80	16/5.7
93CL1525L07	15	25	70	100	16/5.7
93CL1525L09	15	25	90	120	16/5.7
93CL1525L10	15	25	105	135	16/5.7

03 Iliac Extension



Catalog Number	Proximal Ø (mm)	Distal Ø (mm)	Usable Length (mm)	Total Length (mm)	OD Delivery System (F/mm)
93IE1310L05	13	10	50	80	16/5.7
93IE1313L05	13	13	50	80	16/5.7
93IE1913L05	19	13	50	80	16/5.7
93IE1919L05	19	19	50	80	16/5.7
93IE1922L05	19	22	50	80	16/5.7
93IE2219L05	22	19	50	80	16/5.7
93IE2222L05	22	22	50	80	16/5.7
93IE2227L05	22	27	50	80	18/6.4
93IE2722L05	27	22	50	80	18/6.4
93IE2727L05	27	27	50	80	18/6.4

04 Aortic Extension



Catalog Number	Proximal Ø (mm)	Distal Ø (mm)	Covered Length (mm)	OD Delivery System (F/mm)
93AE2323L05	23	23	50	18/6.1
93AE2626L05	26	26	50	18/6.1
93AE2929L05	29	29	50	18/6.1
93AE3232L05	32	32	50	20/6.8
93AE3636L05	36	36	50	20/6.8
93AE3838L05	38	38	50	20/6.8

05 Aorto-Uni-Iliac



93AU2313L11	23	13	105	18/6.1
93AU2613L11	26	13	105	18/6.1
93AU2913L11	29	13	105	18/6.1
93AU3213L11	32	13	105	20/6.8
93AU3613L11	36	13	105	20/6.8

ARTIVION™

Learn more at artivion.com

All products and indications are not available/approved in all markets. All trademarks are owned by Artivion, Inc. or its subsidiaries. On-X Life Technologies, Inc. Jotec GmbH and Ascyrus Medical GmbH are wholly owned subsidiaries of Artivion, Inc. Artivion and E-tegra are Trademarks of Artivion, Inc. or its subsidiaries. © 2022 Artivion, Inc. All rights reserved.

Artivion, Inc.

1655 Roberts Blvd., NW
Kennesaw, GA 30144 USA
Phone: +1-888-427-9654
Fax: +770-590-3753

E-mail: inquiries@artivion.com

For contact information by region, please visit www.artivion.com/contact

JOTEC GmbH, Lotzenäcker 23, 72379 Hechingen, Germany

JT-RS-0930200-EN V01 09/2022